

NOTICE: This document contains references to Varian.
Please note that Varian, Inc. is now part of Agilent
Technologies. For more information, go to
www.agilent.com/chem.



Agilent Technologies

APPLICATION NOTE

Antidiabetics

Application 1996 LC

Column : Polaris 5 μ C18-A, 150 x 4.6mm
Part nr.: A2000150X046
Eluent : MeCN : 75 mM (NH₄)₃PO₄, pH 4.4, 70 : 30
Flow : 1.0 ml/min
Detection : UV @ 254nm

Peak identification

1. glyburide

